



Bilirubin L-Type

Vanadate oxidation method



Ordering information

Product	Package	Code No.
Direct Bilirubin L-Type R1 Buffer Solution	R1: 4 x 70 mL	411-23695
Direct Bilirubin L-Type R2 Vanadate Solution	R2: 4 x 18 mL	413-23895
Total Bilirubin L-Type R1 Buffer Solution	R1: 4 x 70 mL	417-23295
Total Bilirubin L-Type R2 Vanadate Solution	R2: 4 x 18 mL	419-23495
Bilirubin Calibrator	CAL: 4 x for 3 mL	419-73295

Hitachi 7600 / Roche MODULAR

D-type

P-type

Reagent	Total / Direct Bilirubin L-Type	Total / Direct Bilirubin L-Type
ANALYZE		
CH TEST/ TYPE	T-BIL/D-BIL / SERUM	T-BIL/D-BIL / SERUM
ASSAY CODE	2POINT END - 10	2POINT END - 10
ASSAY POINT	22 - 50 - 0 - 0	16 - 34 - 0 - 0
WAVELENGTH(SUB/MAIN)	546 / 450	546 / 450
SAMPLE VOL. (NORMAL)	5.5 - 0.0 - 0	5.5 - 0.0 - 0
(DEC.)		
(INC.)		
DILUENT	H2O - 0	H2O - 0
REAGENT VOL. R1	160 - 0 - () - 0	160 - 0 - () - 0
R2	40 - 0 - () - 0	0 - 0 - () - 0
R3		40 - 0 - () - 0
R4		0 - 0 - () - 0
ABS.LIMIT		
PROZONE LIMIT	0 0 LOWER	0 0 LOWER
CEL DET.	DET. 1	DET. 1
CALIBRATION		
CALIB TYPE	LINEAR	LINEAR
POINT	2 / 2	2 / 2
WEIGHT	0	0
AUTO CALIBRATION	0	0
SD LIMIT	999	999
DUPLICATE LIMIT	99% / 500 Abs	99% / 500 Abs
SENSITIVITY LIMIT	-9999 / 9999	-9999 / 9999
S1 ABS LIMIT	-32000 / 32000	-32000 / 32000
Calibrator	CODE-CONC.-POS.-VOL.-PREDIL	CODE-CONC.-POS.-VOL.-PREDIL
1	*1 - 0.0 - *1 - 5.5 - 0.0	*1 - 0.0 - *1 - 5.5 - 0.0
2	*1 - *1 - *1 - 5.5 - 0.0	*1 - *1 - *1 - 5.5 - 0.0
3		
4		
5		
6		
K-FACTOR		

*1: Input the calibrator code.

Concentration and position of the calibrators are automatically set.

val. Japan 1211D2/kg