



Creatinine M L-Type



METHOD: Enzymatic Creatininase-HMMPS

Ordering Information

Product	Package	Order No.
Creatinine M L-Type R1	R1: 4 x 70 mL	419-08199
Creatinine M L-Type R2	R2: 4 x 24 mL	419-08299

Enzyme Color A (R1) and Enzyme Color B (R2) use as supplied.
 After opening the bottle, this solution is stable for one month at 2 - 10°C.

Measuring Temperature 37°C

AU2700 / AU5400 (BECKMAN COULTER)

Reagent	Creatinine M L-Type
Specific test parameters	
Test Name	CREA*
Type	Serum
Sample: Volume / Dilution	5.0 µL / 0 µL
Reagents: R1 Volume / Dilution	150 µL / 0 µL
R2 Volume / Dilution	50 µL / 0 µL
Wavelength (Pri. / Sec.)	600 / 700
Method	END
Reaction slope	+
Measuring Point 1 (First / Last)	0 / 27
Measuring Point 2 (First / Last)	0 / 10
Linearity	[]%
No-Lag-Time	
Min OD / Max OD (L / H)	
Reagent OD Limit: First (L / H)	-2.0000 / 2.5000
Last (L / H)	-2.0000 / 2.5000
Dynamic range (L / H)	
Correlation Factor (A / B)	1.0000 / 0.0000
Calibration Parameters	
Calibration Type	AB
Formula	$Y = aX + b$
Cuvettes / Counts	
Process	CONC
Calibrator	
Cal. No. - OD - CONC.	
Point 1	* - [] - *
Point 2	* - [] - *
Point 3	* - [] - *
Point 4	* - [] - *
Point 5	* - [] - *
Point 6	* - [] - *
1 -Point Cal. Point	
MB type Factor	

*: To be inserted by the user.

0510D2/vk