FUJIFILM

HA Calibrator Set

Read this instruction sheet thoroughly before use!

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993-71115
1217D2

Intended use

The HA Calibrator Set is designed to be used with the Hyaluronic acid LT reagent for the determination of hyaluronic acid (HA) in serum and plasma.

Reagents and ingredients

(Do not freeze!)

HA Calibrator 1:	(50 ng/mL)	1 bottle x 2 mL	Store at 2–10 °C
HA Calibrator 2:	(100 ng/mL)	1 bottle x 2 mL	Store at 2–10 °C
HA Calibrator 3:	(200 ng/mL)	1 bottle x 2 mL	Store at 2–10 °C
HA Calibrator 4:	(500 ng/mL)	1 bottle x 2 mL	Store at 2–10 °C
HA Calibrator 5:	(1000 ng/mL)	1 bottle x 2 mL	Store at 2–10 °C

Containing sodium hyaluronate from *Streptococcus zooepidemicus* in phosphate buffer.

Reagent preparation

Use the calibrators as supplied.

If the calibrator is stored after opening, tightly cap the bottle and store it at 2–10 °C.

Determination of assigned values

Assigned values: See bottle labels

Determination: The HA concentration was determined by measuring a calibrator with the Hyaluronic acid LT reagent and the Wako 1st HA standard. HA concentration of the Wako 1st HA standard was determined by the carbazole-sulfuric acid method.¹

Warnings and precautions

- For *in vitro* diagnostic use.
- The usage and application of these calibrators are reserved for professional use only. Please refer to respective national and local regulations and legislation.
- Not to be used internally in humans and animals.
- Store the calibrators under the specified conditions. Do not use calibrators past the expiration date stated on each calibrator container label.
- Do not use this calibrator set for any other purpose than described herein. Refer to the package insert for Hyaluronic acid LT.
- Do not use the containers and other materials in the package for any purposes other than those described herein.
- If the calibrators come in contact with mouth, eyes or skin, wash off immediately with a large amount of water. Consult a physician if necessary.
- When discarding the calibrators, dispose of them according to local or national regulations.
- Each calibrator contains 0.09 % sodium azide as a stabilizer. Sodium azide may react with lead or copper plumbing to form explosive compounds. Even though the calibrator contains minute quantities of sodium azide, drains should be flushed well with a large amount of water, when discarding the calibrator.

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