

Product Information

Reagents	Packaging size	Article No.
Hyaluronic acid LT	R1: 2 x 15 mL (HABP Reagent)	992-71185
	R2: 2 x 6 mL (HA Latex Reagent)	
Hyaluronic acid LT	R1: 2 x 31 mL (HABP Reagent)	993-71095
	R2: 2 x 11 mL (HA Latex Reagent)	
HA Calibrator Set	CAL: 5 conc. x 2 mL	993-71115
HA Control Set	CONTROL: 2 x 2 conc. x 2 mL	998-71165

Use R1 and R2 as supplied. After opening bottle, store at 2 - 10°C and use within one month.

Hitachi 912

Reagent	Hyaluronic acid LT
Analyze	
Test	HA
Assay / Time	2 Point End - 10
Assay Points	19 - 28 - 0 - 0
Wavelength (Sub / Main)	None / 800
Dilution	00001 / 99
Sample Volume (Normal)	4 - 0 - 0
(Decrease), (Increase)	
Abs. Limit	32000 / Increase
Prozone Limit	-32000 / Lower
Reagent R1	240
R2	
R3	80
R4	
Calibration	
Calibration Type	Spline
Number of Points	6
Span Point	5
Calib. Weight	0
Auto Calibration	
SD Limit	999.9
Duplicate Limit	99% 32000 Abs.
Sensitivity Limit	-99999 - 999999
S1 Ext Limit	-32000 - 32000
Range	
Data Mode	Inactiv
Unit	ng / mL
Instrument Factor	(Y = aX + b)
a	
b	
Compensated	-99999 - 999999
Standards	
	Conc. - Pos. - Vol. - Dil
1	(0.0) *1 - *2 - 4 - 0
2	*1 - *2 - 4 - 0
3	*1 - *2 - 4 - 0
4	*1 - *2 - 4 - 0
5	*1 - *2 - 4 - 0
Calibration List	
K-Factor	
System Parameter	
Linearity Check (4 - 8)	
(9 -)	

*1: Assigned value of the calibrator

*2: Position of the calibrator

0414D3/vk