

Ordering Information

Reagents	Package	Order No.
HDL-C L-Type R1	R1: 4 x 60 mL	412-72395
HDL-C L-Type R2	R2: 4 x 20 mL	412-72495
HDL-C L-Type R1	R1: 4 x 270 mL	418-72395
HDL-C L-Type R2	R2: 4 x 90 mL	418-72495
HDL-C Calibrator	CAL: 4 x for 3 mL	416-51095

Hitachi 7600/ Roche MODULAR

D-type

P-type

Reagent	HDL-C L-Type	HDL-C L-Type
ANALYZE		
CH TEST/ TYPE	HDL / SERUM	HDL / SERUM
ASSAY	2POINT END - 10	2POINT END - 10
POINT	16 - 34 - 0 - 0	16 - 34 - 0 - 0
WAVELENGTH(SUB/MAIN)	700 / 600	700 / 600
SAMPLE VOL. (NORMAL)	3.0 - 0.0 - 0	3.0 - 0.0 - 0
(DEC.)	2.0 - 0.0 - 0	2.0 - 0.0 - 0
(INC.)	3.0 - 0.0 - 0	3.0 - 0.0 - 0
DILUENT	H2O - 0	H2O - 0
REAGENT VOL. R1	200 - 0 - () - 0	200 - 0 - () - 0
R2	67 - 0 - () - 0	0 - 0 - () - 0
R3		67 - 0 - () - 0
R4		0 - 0 - () - 0
ABS LIMIT	INCREASE	INCREASE
PROZONE LIMIT	0 0 LOWER	0 0 LOWER
CELL DET.	DET. 1	DET. 1
CALIBRATION		
CALIB. TYPE	LINEAR	LINEAR
POINT	2 / 2	2 / 2
WEIGHT	0	0
UNIT	mg/dL	mg/dL
AUTO CALIBRATION	0	0
SD LIMIT	0.1	0.1
DUPLICATE LIMIT	5% / 80 Abs	5% / 80 Abs
SENSITIVITY LIMIT	10 / 35	10 / 35
S1 ABS LIMIT	-32000 / 32000	-32000 / 32000
CALIBRATOR	Code CONC.POS.VOL.PREDIL	Code CONC.POS.VOL.PREDIL
1	*1 *1(0.0) *1 3.0 0.0	*1 *1(0.0) *1 3.0 0.0
2	*1 *1 *1 3.0 0.0	*1 *1 *1 3.0 0.0
3		
4		
5		
6		
K-FACTOR		

*1: Input the Calibrator Code

Concentration and position of calibrators are automatically set