

Ordering information

Product	Package	Code No.
Direct Bilirubin L-Type R1 (Buffer Solution)	R1: 4 x 70 mL	411-23695
Direct Bilirubin L-Type R2 (Vanadate Solution)	R2: 4 x 18 mL	413-23895
Total Bilirubin L-Type R1 (Buffer Solution)	R1: 4 x 70 mL	417-23295
Total Bilirubin L-Type R2 (Vanadate Solution)	R2: 4 x 18 mL	419-23495
Bilirubin Calibrator	CAL: 4 x for 3 mL	419-73295

Beckman Coulter AU2700/ AU5400

Reagent	Bilirubin L-Type
Specific test parameters	Page 1 / 2
Test No. / Name	T-BIL / D-BIL
Sample	Serum
Sample Vol / Dil	5.5 / 0
Reagent 1 Vol / Dil	160 / 0
Reagent 2 Vol / Dil	40 / 0
Wave (Main / Sub)	450 / 540
Method	END
Reaction	-
Point 1 (1st / Lst)	0 / 27
Point 2 (1st / Lst)	0 / 10
Linearity (1st / 2nd)	
No-Lag-Time	1
Min OD / Max OD (L/H)	
Reagent OD (1stL / 1stH)	-2.0000 / 2.5000
Limit (LstL / LstH)	-2.0000 / 2.5000
Dynamic range (L / H)	
Correction Factor (A / B)	1.0000 / 0.0000
Calibration Parameters	
Calibration Type	AB
Formula	Y=AX+B
Cuvettes	
process	conc.
Calibrator	Cal # - OD - Conc- Factor
Point 1	*1 - [] - 0 - 9999999
Point 2	*1 - [] - *2 - 9999999

*1: To be inserted by the user

*2: Use calibrator concentration