



# CRP-HSII-LT



METHOD: Turbidimetric Immunoassay

## Ordering information

Reagent	Package	Code No.
CRP-HSII-LT (Buffer)	R1: 2 x 40 mL	462-39493
CRP-HSII-LT (Latex Reagent)	R2: 2 x 40 mL	468-39593
CRP-HSII-LT Calibrator Set HO	CAL: 5 conc. x 2 mL	460-25801

Measuring Temperature 37°C

## AU640 (BECKMAN COULTER)

Reagent	CRP-HSII-LT
<b>Specific test parameters</b>	
Test No. / Name	* / CRP
Sample	Serum
Operation	YES
Sample Vol / Dil	3.5 / 0
Reagent 1 Vol / Dil	150 / 0
Reagent 2 Vol / Dil	100 / 0
Wave (Main / Sub)	600 / -
Method	END
Reaction	+
Point 1 (1st / Lst)	0 / 26
Point 2 (1st / Lst)	0 / 11
Linearity (1st / 2nd)	
No-Lag-Time	
Min OD / Max OD (L / H)	
Sample Dilution	1
Reagent OD (1stL / 1stH)	-2.0000 / 2.5000
Limit (LstL / LstH)	-2.0000 / 2.5000
Dynamic range (L / H)	
Correction Factor (A / B)	1.0000 / 0.0000
<b>Calibration Parameters</b>	
Calibration Type	6AB
Formula	SPLINE
Cuvettes	
Process	CONC
<b>Calibrator</b>	<b>Cal # - OD - Conc- Factor</b>
Point 1	* - [ ] - * -
Point 2	* - [ ] - * -
Point 3	* - [ ] - * -
Point 4	* - [ ] - * -
Point 5	* - [ ] - * -
Point 6	* - [ ] - * -
Point 7	
1 point cal point	
MB type Factor	
<b>Data Check Condition</b>	
<b>Data Check-3</b>	
Check point 1	12
Step interval	3
Value 1	2,0000
Value 2	0,02
Fast Point	12
Last Point	14

\*1: To be inserted by the user.

0111D1/vk

\*2: Assigned value of the calibrator (mg/L or mg/dL).