

## Ordering information

Code No.	Product	Package
415-73301	Autokit Total Ketone Bodies R1 Set (Autokit T-KB R1 Set)	R1a: 2 x 27 mL R1b: 2 x for 27 mL
413-73601	Autokit 3-HB R2 Set	R2a: 2 x 9 mL R2b: 2 x for 9 mL
412-73791	Ketone Body Calibrator • 300 (3-HB: 300µmol/L)	CAL: 4 x 5 mL

## Hitachi 7070 (7060) / 911 (904) / 912

[High sensitivity method]

TEST	T-KB (Total ketone bodies)	T-KB (Total ketone bodies)
ASSAY CODE	RATE-A : 10	RATE-A : 10
WAVELENGTH (SUB/MAIN)	600 / 415	600 / 415
ASSAY POINT	19 - 23 - 0 - 0	19 - 23 - 0 - 0
DILUTION	/ 0	/ 0
SAMPLE VOL. (NORM.)	4 - 0 - 0	17 - 0 - 0
(DEC.)		
(INC.)		
ABS.LIMIT	13000 / INC.	13000 / INC.
PROZONE LIMIT	-32000 / LOWER	-32000 / LOWER
REAGENT R1	270 - 0 - ( ) - 0	270 - 0 - ( ) - 0
R2	0 - 0 - ( ) - 0	0 - 0 - ( ) - 0
R3	90 - 0 - ( ) - 0	90 - 0 - ( ) - 0
R4	0 - 0 - ( ) - 0	0 - 0 - ( ) - 0
CALIB.TYPE	LINEAR - 2 - 2 - 0	LINEAR - 2 - 2 - 0
AUTO CALIB.	0	0
SD LIMIT		
DUPLICATE LIMIT	1000	1000
SENSITIVITY LIMIT	0	0
S1 ABS LIMIT	-32000 / 32000	-32000 / 32000
COMPENSATED LIMIT		
UNIT	mmol/L	mmol/L
INSTRUMENT a	1.0	1.0
FACTOR (Y=aX+b) b	0.0	0.0
<b>STANDARD</b>	<b>Conc.Pos.Vol.Dil</b>	<b>Conc.Pos.Vol.Dil</b>
1	0.0 - W - 4 - 0 - 0	0.0 - W - 17 - 0 - 0
2	*1 - *2 - 4 - 0 - 0	*1 - *2 - 17 - 0 - 0
3		
4		
5		
6		
<b>CALIBRATION LIST</b>		
K-FACTOR		
<b>SYSTEM PARAMETER</b>		
LINEARITY CHECK (4 - 8)		
(9 - )		

\*1: Assigned value of a calibrator

\*2: Position of the calibrator