

FUJIFILM

CRP-HS LT Calibrator Set

Read this instruction sheet thoroughly before use!



460-25801
0318D2

Intended use

The CRP-HS LT Calibrator Set is designed to be used with the CRP-HS II LT reagent for the determination of C-Reactive Protein (CRP).

Reagents and ingredients

There are two different Calibrator Sets available:

CRP-HS LT Calibrator Set HO 5 concentration x 2 mL
(for automated analyzers like Hitachi, AU and Advia Systems)

Product	Concentration	Package
CRP-HS LT Calibrator HO 1	0.2 mg/dL	CAL 1: 1 x 2 mL
CRP-HS LT Calibrator HO 2	1.0 mg/dL	CAL 2: 1 x 2 mL
CRP-HS LT Calibrator HO 3	4.0 mg/dL	CAL 3: 1 x 2 mL
CRP-HS LT Calibrator HO 4	18.0 mg/dL	CAL 4: 1 x 2 mL
CRP-HS LT Calibrator HO 5	35.0 mg/dL	CAL 5: 1 x 2 mL

CRP-HS LT Calibrator Set T 7 concentration x 2 mL
(for automated analyzers like Architect and Aeroset Systems)

Product	Concentration	Package
CRP-HS LT Calibrator T 1	0.3 mg/dL	CAL 1: 1 x 2 mL
CRP-HS LT Calibrator T 2	0.8 mg/dL	CAL 2: 1 x 2 mL
CRP-HS LT Calibrator T 3	3.5 mg/dL	CAL 3: 1 x 2 mL
CRP-HS LT Calibrator T 4	8.0 mg/dL	CAL 4: 1 x 2 mL
CRP-HS LT Calibrator T 5	16.0 mg/dL	CAL 5: 1 x 2 mL
CRP-HS LT Calibrator T 6	28.0 mg/dL	CAL 6: 1 x 2 mL
CRP-HS LT Calibrator T 7	35.0 mg/dL	CAL 7: 1 x 2 mL

Determination of assigned values

The CRP concentration was determined by Latex Immunoassay using ERM-DA 470 (IRMM) as primary standard.

Reagent preparation

Use the calibrators as supplied. After opening the calibrators, it is recommended to use them immediately. When the opened calibrators are stored, cap the bottles to keep from drying out and store them at 2–10 °C.

Warnings and precautions

- For *in-vitro*-diagnostic use.
- The usage and application of this test is reserved for professional use only.
- Not to be used internally in humans or animals.
- Do not use the calibrators for any other procedures than described herein. Refer to the package insert for CRP-HS II LT reagent before use of the calibrators!
- Do not use the calibrators for any other purpose than described herein.
- Do not use the containers and other materials in the package for any purposes other than those described herein.
- Store the calibrators under the specified conditions. Do not use calibrators past the expiration date stated on each calibrator container label.
- If the calibrators come in contact with mouth, eyes or skin, wash off immediately with a large amount of water. Consult a physician if necessary.
- The calibrators are prepared from humansera which were found to be negative for HbsAg, anti-HIV antibody and anti-HCV antibody. Because no test can offer complete assurance that products derived from human blood will not transmit infectious agents, it is recommended that the calibrators should be handled with the same precautions used for patient specimens.
- The calibrators contain 0.09 % sodium azide as a stabilizer. Sodium azide may react with copper or lead plumbing to form explosive compounds. Even though the reagent contains minute quantity of sodium azide, drains should be flushed well with a large amount of water, when discarding the reagent.
- When discarding the calibrators, dispose of them according to local or national regulations

CRP-HS LT Calibrator Set contains components classified as follows according to Regulation (EC) No. 1272/2008: 5-chloro-2-methyl-2H-isothiazol-3-one [EC No. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC No. 220-239-6] (3:1).

Hazard pictogram



Signal word Warning

Hazard statements

H317 May cause an allergic skin reaction.

Precautionary statements

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
 P262 Do not get in eyes, on skin, or on clothing.
 P281 Use personal protective equipment as required.
 P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
 P302+P352 IF ON SKIN: Wash with plenty of soap and water.

Storage condition and shelf life

(1) **Storage condition:** 2–10 °C.
 (2) **Shelf life:** 1 year after production.

Quality control

A quality control program is recommended for all clinical laboratories.

Ordering information

Code No.	Product	Package
460-25801	CRP-HS LT Calibrator Set HO (for automated analyzers like Hitachi, AU and Advia Systems)	CAL: 5 conc. x 2 mL
on request	CRP-HS LT Calibrator Set T (for automated analyzer like Architect and Aeroset Systems)	CAL: 7 conc. x 2 mL

Manufacturer:

FUJIFILM Wako Chemicals Europe GmbH
 Fuggerstr. 12 · 41468 Neuss · Germany
 Telephone: +49 2131 311 272 · Fax: +49 2131 311 110
 diagnostics_wkeu@fujifilm.com · www.wako-chemicals.de

Wako