

Ordering Information

| Reagents | Package | Order No. |
|---|-------------------|-----------|
| HDL-C L-Type / Immuno-inhibition Method | | |
| HDL-C L-Type R1 | R1: 4 x 60 mL | 412-72395 |
| HDL-C L-Type R2 | R2: 4 x 20 mL | 412-72495 |
| HDL-C L-Type R1 | R1: 4 x 270 mL | 418-72395 |
| HDL-C L-Type R2 | R2: 4 x 90 mL | 418-72495 |
| HDL-C Calibrator | CAL: 4 x for 3 mL | 416-51095 |
| LDL-C L-Type / Enzymatic Selective Protection Method | | |
| LDL-C L-Type R1 | R1: 4 x 60 mL | 419-24017 |
| LDL-C L-Type R2 | R2: 4 x 20 mL | 419-24027 |
| LDL-C L-Type R1 | R1: 1 x 500 mL | 991-24015 |
| LDL-C L-Type R2 | R2: 1 x 500 mL | 991-24025 |
| LDL-C Calibrator | CAL: 4 x for 1 mL | 419-24032 |

AU5400 (Beckman Coulter)

| Reagent | HDL-C L-Type | LDL-C L-Type |
|---------------------------------|----------------------------------|----------------------------------|
| Specific test parameters | | |
| Test No. / Name | * / HDL-C | * / LDL-C |
| Sample | Serum | Serum |
| Sample Vol / Dil / ratio | 1.6 / 0 / 1 | 1.6 / 0 / 1 |
| Reagent 1 Vol / Dil | 100 / 35 | 120 / 10 |
| Reagent 2 Vol / Dil | 34 / 0 | 40 / 0 |
| Wave (Main / Sub) | 600 / 700 | 600 / 700 |
| Method | END | END |
| Reaction | + | + |
| Point 1 (1st / Lst) | 10 / 27 | 0 / 27 |
| Point 2 (1st / Lst) | | 0 / 10 |
| Linearity % | | |
| No-Lag-Time | | |
| Min OD / Max OD (L / H) | | |
| Reagent OD (1stL/1stH) | -2.0000 / 2.5000 | -2.0000 / 2.5000 |
| (LstL/LstH) | -2.0000 / 2.5000 | -2.0000 / 2.5000 |
| Dynamic range (L / H) | | |
| Correction Factor (A / B) | 1.000000 / 0.000000 | 1.000000 / 0.000000 |
| Calibration Parameters | | |
| Calibration Type | AB | AB |
| Formula | 1 (Y=AX+B) | 1 (Y=AX+B) |
| Seq. No | 2 | 2 |
| Process | | |
| Calibrator | Cal # - OD - Conc- Factor | Cal # - OD - Conc- Factor |
| Point 1 | * - [] - * - | * - [] - * - |
| Point 2 | | |
| Point 3 | | |
| Point 4 | | |
| Point 5 | | |
| Point 6 | | |
| Point 7 | | |

* : To be inserted by the user