# FUJIFILM

## Bilirubin Calibrator

Read this instruction sheet thoroughly before use!

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419-73295  
1223D6

### Intended use

Wako Bilirubin Calibrator is designed to be used for the calibration with Wako's Total and Direct bilirubin reagents of the vanadate oxidation method in the determination of total and direct bilirubin.

### Reagents

Bilirubin CAL: 4 bottles x for 3 mL.

Contains bilirubin from bovine bile, ditaurobilirubin, and 0.005% (when reconstituted) sodium azide. Store at 2–10 °C. Protect from light.

### Assigned values and their determination

Assigned values: See bottle labels.

Determination: Based on the US NIST bilirubin standard material (SRM916) with the vanadate oxidation method.

### Reagent preparation

Accurately add 3.0 mL of distilled or deionized water to dissolve the contents of each of the Calibrator. The reconstituted Calibrator is stable for 5 days at 2 – 10 °C under protection from light.

### Warnings and precautions

- For *in vitro* diagnostic use.
- Not to be used internally in humans or animals.
- Do not use this calibrator described above in any procedures other than those described herein.
- Do not use reagents past the expiration date stated on each reagent container label.
- The assigned values on the label are valid for only the vanadate oxidation method, and can not be used with other methods because the reactivity of the other methods to ditaurobilirubin is different from that of the vanadate method.
- Store the reconstituted solution in a refrigerator or a freezer under protection from light. Do not repeat freeze and thaw more than once.
- This product contains 0.005 % (when reconstituted) sodium azide as a stabilizer. Sodium azide may react with copper or lead plumbing to form explosive compounds. Even though this product contains minute quantities of sodium azide, drains should be flushed well with a large amount of water when discarding the solution.
- Be careful not to cut yourself with the aluminium cap when you remove it from the vial.
- The vial is stoppered at reduced pressure. Slowly remove the stopper in order not to blow the powder in the vial.
- If this calibrator comes in contact with mouth, eye or skin, wash off immediately with a large amount of water. Consult a physician if necessary.
- When discarding this calibrator, dispose of it according to local or national regulations.

### Labelling according to Regulation (EC) No 1272/2008

- The product is classified and labelled according to the CLP regulation.
- The product does not contain a SVHC according to REACH, Article 57 >0,1 %

### Hazard-determining components of labelling

- Globulins,  $\gamma$ -

### Hazard pictograms

Warning

### Hazard statements

- Causes serious eye irritation.
- May cause respiratory irritation.

### Precautionary statements

- Avoid breathing dust/fume/gas/mist/vapours/spray.
- Wash thoroughly after handling.
- Wear eye protection/face protection.
- IF INHALED: Remove person to fresh air and keep comfortable for breathing.
- IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- Call a POISON CENTER/doctor if you feel unwell.
- Contains Gentamicin, sulfate (salt). May produce an allergic reaction.

### Ordering information

Code No.	Product	Package
419-73295	Bilirubin Calibrator	CAL: 4 x for 3 mL

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**FUJIFILM Wako Chemicals Europe GmbH**