#### Characteristics

Control Serum I/Control Serum II is a control serum for the normal/abnormal range. Control Serum I/Control Serum II is a lyophilised serum based on human serum suitable for use with both manual and automated analytical procedures to monitor the precision of laboratory procedures.

#### Summary and principle

The use of quality control materials is indicated as an objective assessment of the precision of methods and techniques in use and is an integral part of good laboratory practices.

#### Reagent

This product is prepared from human serum with added constituents of human and animal origin.

#### Actual values and range

The mean values printed in this insert were derived from replicate analyses and are specific for this lot. The tests listed were performed by the reagent manufacturer. Individual laboratory means should fall within the corresponding acceptable range, however, laboratory means may vary from the listed values during the life of this product. It is recommended that each laboratory establishes its own means and acceptable ranges and uses those provided only as guides. See data sheet.

#### Preparation

Add exactly 5 mL dist. water and let stand for 30 min. Dissolve the contents completely by swirling, avoiding the formation of foam. Do not shake. Before sampling, allow the control to reach room temperature. This product should be treated the same as patient specimens and run in accordance with the instructions accompanying the instrument, kit or reagent being used.

#### Storage and Stability

This product will be stable until the expiration date when stored unopened at 2 to 8 °C. Once the control is opened the analytes will be stable for one day at 15–20 °C, 7 days at 2–8 °C, and one month at -20 °C (see remarks).

#### Warnings and precautions

- This product should not be used past the expiration date.
- If there is evidence of microbiological contamination or excessive turbidity in the product, discard the vial.
- Dispose of any discarded material in accordance with the requirements of your local waste management authorities.
- This product is not intended for use as a calibrator.

#### Remarks

- Acid phosphatase has to be measured immediately after dissolving.
- Stability of the CPK: at 2–20 °C, 8 hours.
- Dissolved serum has to be kept in the dark. CK and bilirubin activities decrease by light action.
- Shorter stability characteristics may apply to some enzymes, phosphorus, triglycerides and bilirubin.
- Wako Control Serum I/Control Serum II has been prepared from normal human sera that were found non-reactive against HIV-1, 2 antibodies and hepatitis C antibodies as well as HBs antigen using FDA approved methods. Since the risk of infection cannot be excluded with certainty, however, the products must be handled as potentially biohazardous like patients' sera.

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