Read this instruction sheet thoroughly before use!

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998-71165
1217D2

Intended use

The HA Control Set is designed to be used as a quality control material for the quantitative determination of hyaluronic acid (HA) using the Hyaluronic acid LT reagent.

Reagents and ingredients

HA Control 1: 2 bottles x 2 mL Store at 2–10 °C (Do not freeze)
HA Control 2: 2 bottles x 2 mL Store at 2–10 °C (Do not freeze)

Containing sodium hyaluronate from *Streptococcus zooepidemicus* in phosphate buffer.

Reagent preparation

Use the controls as supplied.

If the control is stored after opening, tightly cap the bottle and store it at 2–10 °C.

Determination of assigned values

Assigned values: See bottle labels

Determination: The HA concentration was determined by measuring a control with the Hyaluronic acid LT reagent and HA calibrators whose HA levels are traceable to the Wako 1st HA standard. HA concentration of the Wako 1st HA standard was determined by the carbazole-sulfuric acid method.¹

Warnings and precautions

- For *in vitro* diagnostic use.
- The usage and application of these controls is reserved for professional use only. Please refer to respective national and local regulations and legislation.
- Not to be used internally in humans and animals.
- Store the controls under the specified conditions. Do not use controls past the expiration date stated on each control container label.
- Do not use this control set for any other purpose than described herein. Refer to the package insert for Hyaluronic acid LT.
- Do not use the containers and other materials in the package for any purposes other than those described herein.
- If the controls come in contact with mouth, eyes or skin, wash off immediately with a large amount of water. Consult a physician if necessary.
- When discarding the controls, dispose of them according to local or national regulations.
- Each control contains 0.09 % sodium azide as a stabilizer. Sodium azide may react with lead or copper plumbing to form explosive compounds. Even though the control contains minute quantities of sodium azide, drains should be flushed well with a large amount of water, when discarding the control.

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