

Ordering information

Reagents	Package	Order No.
LDL-C L-Type R1	R1: 4 x 60 mL	419-24017
LDL-C L-Type R2	R2: 4 x 20 mL	419-24027
LDL-C L-Type R1	R1: 1 x 500 mL	991-24015
LDL-C L-Type R2	R2: 1 x 500 mL	991-24025
LDL-C Calibrator	CAL: 4 x for 1 mL	419-24032

AU400 / AU600 / AU640 (BECKMAN COULTER)

Reagent	LDL-C L-Type
Specific Test Parameters	
Test Name	*1 / LDL-C
Type	Serum
Sample: Volume / Dilution	2.5 µL / 0 µL
Reagents: R1 Volume / Dilution	210 µL / 0 µL
R2 Volume / Dilution	70 µL / 0 µL
Wavelength (Pri. / Sec.)	600 / 700
Method	END
Reaction slope	+
Measuring Point 1 (First / Last)	0 / 27
Measuring Point 2 (First / Last)	0 / 10
Linearity	[]%
No-Lag-Time	
Min OD / Max OD (L / H)	
Reagent OD Limit: First (L / H)	-2.0000 / 2.5000
Last (L / H)	-2.0000 / 2.5000
Dynamic Range (L / H)	0.02 - 10.30
Correlation Factor (A / B)	1.000000 / 0.000000
Calibration Specific	
Calibration Type	AB
Formula	1 (Y = aX +b)
Cuvette / Counts	2 / *1
Process	CONC
Calibrator	CAL# - OD - CONC
Point 1	*1 - [] - *2
Point 2	
Point 3	
Point 4	
Point 5	
Point 6	
Point 7	
1-Point Cal. Point	
MB Type Factor	

*1: To be inserted by the user

*2: Input the concentration of the calibrator