

### Ordering Information

Reagents	Package	Order No.
LDL-C L-Type R1	R1: 4 x 60 mL	419-24017
LDL-C L-Type R2	R2: 4 x 20 mL	419-24027
LDL-C L-Type R1	R1: 1 x 500 mL	991-24015
LDL-C L-Type R2	R2: 1 x 500 mL	991-24025
LDL-C Calibrator	CAL: 4 x for 1 mL	419-24032

### Hitachi 7600/ Roche MODULAR

#### D-type

#### P-type

Reagent	LDL-C L-Type	LDL-C L-Type
<b>ANALYZE</b>		
CH TEST/ TYPE	LDL / SERUM	LDL / SERUM
ASSAY	2POINT END - 10	2POINT END - 10
POINT	16 - 34 - 0 - 0	16 - 34 - 0 - 0
WAVELENGTH(SUB/MAIN)	700 / 600	700 / 600
SAMPLE VOL. (NORMAL)	2.0 - 0.0 - 0	2.0 - 0.0 - 0
(DEC.)	2.0 - 0.0 - 0	2.0 - 0.0 - 0
(INC.)	2.0 - 0.0 - 0	2.0 - 0.0 - 0
DILUENT	H2O - 0	H2O - 0
REAGENT VOL. R1	180 - 0 - ( ) - 0	180 - 0 - ( ) - 0
R2	60 - 0 - ( ) - 0	0 - 0 - ( ) - 0
R3		60 - 0 - ( ) - 0
R4		0 - 0 - ( ) - 0
ABS LIMIT	INCREASE	INCREASE
PROZONE LIMIT	0 0 LOWER	0 0 LOWER
CELL DET.	DET. 1	DET. 1
<b>CALIBRATION</b>		
CALIB. TYPE	LINEAR	LINEAR
POINT	2 / 2	2 / 2
WEIGHT	0	0
UNIT	mg/dL	mg/dL
AUTO CALIBRATION	0	0
SD LIMIT	999	999
DUPLICATE LIMIT	99% / 500 Abs	99% / 500 Abs
SENSITIVITY LIMIT	-9999 / 9999	-9999 / 9999
S1 ABS LIMIT	-32000 / 32000	-32000 / 32000
<b>CALIBRATOR</b>	Code CONC.POS.VOL.PREDIL	Code CONC.POS.VOL.PREDIL
1	*1 *1(0.0) *1 2.0 0.0	*1 *1(0.0) *1 2.0 0.0
2	*1 *1 *1 2.0 0.0	*1 *1 *1 2.0 0.0
3		
4		
5		
6		
K-FACTOR		

\*1: Input the Calibrator Code.

Concentration and position of calibrators are automatically set