

Ordering information

| Product | Package | Code No. |
|---|-------------------|-----------|
| Direct Bilirubin L-Type R1 Buffer Solution | R1: 4 x 70 mL | 411-23695 |
| Direct Bilirubin L-Type R2 Vanadate Solution | R2: 4 x 18 mL | 413-23895 |
| Total Bilirubin L-Type R1 Buffer Solution | R1: 4 x 70 mL | 417-23295 |
| Total Bilirubin L-Type R2 Vanadate Solution | R2: 4 x 18 mL | 419-23495 |
| Bilirubin Calibrator | CAL: 4 x for 3 mL | 419-73295 |

Hitachi 7600 / Roche MODULAR

D-type

P-type

| Reagent | Total / Direct Bilirubin L-Type | Total / Direct Bilirubin L-Type |
|----------------------|------------------------------------|------------------------------------|
| ANALYZE | | |
| CH TEST/ TYPE | T-BIL / D-BIL / SERUM | T-BIL / D-BIL / SERUM |
| ASSAY CODE | 2POINT END - 10 | 2POINT END - 10 |
| ASSAY POINT | 22 - 50 - 0 - 0 | 16 - 34 - 0 - 0 |
| WAVELENGTH(SUB/MAIN) | 546 / 450 | 546 / 450 |
| SAMPLE VOL. (NORMAL) | 5.5 - 0.0 - 0 | 5.5 - 0.0 - 0 |
| (DEC.) | | |
| (INC.) | | |
| DILUENT | H2O - 0 | H2O - 0 |
| REAGENT VOL. R1 | 160 - 0 - () - 0 | 160 - 0 - () - 0 |
| R2 | 40 - 0 - () - 0 | 0 - 0 - () - 0 |
| R3 | | 40 - 0 - () - 0 |
| R4 | | 0 - 0 - () - 0 |
| ABS.LIMIT | | |
| PROZONE LIMIT | 0 0 LOWER | 0 0 LOWER |
| CEL DET. | DET. 1 | DET. 1 |
| CALIBRATION | | |
| CALIB TYPE | LINEAR | LINEAR |
| POINT | 2 / 2 | 2 / 2 |
| WEIGHT | 0 | 0 |
| AUTO CALIBRATION | 0 | 0 |
| SD LIMIT | 999 | 999 |
| DUPLICATE LIMIT | 99% / 500 Abs | 99% / 500 Abs |
| SENSITIVITY LIMIT | -9999 / 9999 | -9999 / 9999 |
| S1 ABS LIMIT | -32000 / 32000 | -32000 / 32000 |
| Calibrator | CODE-CONC.-POS.-VOL.-PREDIL | CODE-CONC.-POS.-VOL.-PREDIL |
| 1 | *1 - 0.0 - *1 - 5.5 - 0.0 | *1 - 0.0 - *1 - 5.5 - 0.0 |
| 2 | *1 - *1 - *1 - 5.5 - 0.0 | *1 - *1 - *1 - 5.5 - 0.0 |
| 3 | | |
| 4 | | |
| 5 | | |
| 6 | | |
| K-FACTOR | | |

*1: Input the calibrator code
Concentration and position of the calibrators are automatically set