



Apolipoprotein A1-HA



Automated turbidimetric immunoassay

Ordering information

Order No.	410-27201
Product	Apolipoprotein A1-HA
Package	R1: 3 x 45 mL R2: 1 x 19 mL

Hitachi 902 / 904 / 911 / 912 / Assay Parameters

REAGENT	APO A1-HA
TEST	APO A1
ASSAY CODE	2POINT END - 10
WAVELENGTH (SUB)	0
WAVELENGTH (MAIN)	700
ASSAY POINT	15 - 31 - 0 - 0
DILUTION	/ 0
SAMPLE VOL. (NORM.)	3 - 0 - 0
(DEC.)	
(INC.)	
ABS.LIMIT	
PROZONE LIMIT	-32000 LOWER
REAGENT R1	350 - 0 - () - 0
R2	0 - 0 - () - 0
R3	50 - 0 - () - 0
R4	0 - 0 - () - 0
CALIB.TYPE	LINEAR - 2 - 2 - 0
AUTO CALIB.	
SD LIMIT	999.9
DUPLICATE LIMIT	500
SENSITIVITY LIMIT	0
S1 ABS LIMIT	-32000/32000
UNIT	mg/dL
INSRUMENT FACTOR	(Y = aX + b)
a	1.0
b	0
STANDARD	Conc. - Pos. - Vol. - Dil.
1	0 - 1 - 3 - 0 - 0
2	*1 - *2 - 3 - 0 - 0
3	
4	
5	
6	

*1: Concentration of standard

0506D1/vk

*2: Position of standard